




Note: The most current version of these instructions can be found online at www.grahamfield.com

- ⚠ WARNING: Important! Read and understand these instructions before assembling or using the Sit-to-Stand Buttock Strap. If you do not understand any part of these warnings, cautions or instructions, contact a healthcare professional for direction in the use of this product. If the Sit-to-Stand Buttock Strap is not properly applied, personal injury and damage to the Sit-to-Stand Buttock Strap could result.**
- ⚠ WARNING: If components are damaged or missing, contact your dealer immediately.**
- ⚠ WARNING: Carefully inspect the Buttock Strap weekly for wear and damage to seams, fabric, straps, and strap loops. Torn, cut, frayed or broken Buttock Straps can fail, resulting in serious personal injury to the user. Use only Buttock Straps that are in good condition. Discard and destroy old, unusable Buttock Straps.**
- ⚠ WARNING: The Sit-to-Stand Buttock has a maximum weight capacity of 400 lb. DO NOT exceed the maximum weight capacity.**
- ⚠ WARNING: GF Health Products, Inc. assumes no responsibility for any damage or injury caused by improper application or use of this product.**

APPLICATION

1. Check the patient's weight. Ensure that the patient's weight does not exceed the Buttock Strap's maximum weight capacity.
 2. Review the following procedures carefully before attempting lift; call Graham-Field Technical Support at 800-347-5678 with any questions.
- ⚠ WARNING: The Buttock Strap must be used with the Sit-to-Stand Sling, as shown in Figure 4. The Sit-to-Stand Sling and Buttock Strap are designed to complement each other in providing the greatest comfort and security for patient transfer.**
3. Use long loops, center loops, or short loops as shown in Figure 2, to position Buttock Strap to provide the greatest support to the patient's middle and lower body during the following steps.
 4. Attach Buttock Strap strap **A** to Sit-to-Stand Lift mast lower loop **a** as shown in Figures 1 and 3.
 5. Bring Buttock Strap around the back of the patient's buttocks and attach Buttock Strap strap **B** to Sit-to-Stand Lift mast lower loop **b** as shown in Figures 1 and 3.
 6. Do not lock caster brakes, in order to allow the lift to properly balance weight.
 7. Proceed with lift and transfer as shown in Figure 4.



MAINTENANCE

To clean Buttock Strap, use a mild soap and water solution and a clean cloth. Rinse with clear water, then allow to air dry.

▲ **Caution:** Never use solvents, abrasive cleaners, or scouring pads on any part of the Buttock Strap.

WARRANTY


GF Health Products, Inc. offers a one-year limited warranty against manufacturer's defects on the Sit-to-Stand Buttock Strap. If a product is deemed to be under warranty, GF Health Products, Inc. shall provide, at its option, (1) replacement of any defective part or product or (2) a credit of the original selling price made to GF Health Products, Inc.'s initial customer. The warranty does not include any labor charges incurred in replacement part(s) installation or any associated freight or shipping charges to GF Health Products, Inc.

FEATURES

Sit-to-Stand Buttock Strap

Designed for use with
LF2020 STS Lift. 400 lb
maximum weight capacity.

DSLSA9 1 ea



ITEM	Torso Girth	Body Length	Total Width	Total Length
DSLSA9	34 1/4"	7 1/2"	71 1/4"	5 1/2"



GRAHAM-FIELD.

GF Health Products, Inc. is not responsible for typographical errors. Packaging, warranties, products and specifications are subject to change without notice.

Graham-Field and Lumex are registered trademarks of GF Health Products, Inc.

© October 2009 GF Health Products, Inc. • 2935 Northeast Parkway • Atlanta, GA 30360 • tel 800-347-5678 • fax 800-721-0601 • www.grahamfield.com