

Connect with a NewLeaf Home Medical Advisor at 888.739.4330 for personal assistance
or to order your very own Invacare micorAIR mattress today!



Shop at [NewLeafHomeMedical.com](https://www.NewLeafHomeMedical.com)

INVACARE® MICROAIR® MA500
THERAPEUTIC SUPPORT SURFACE



microAIR®

Invacare® microAIR® MA500

Alternating Pressure Low Air Loss Mattress System

Low - Medium Risk Stage I-III

A lightweight and easy-to-use Alternating Pressure Mattress System with Low Air Loss aids in the Prevention and Treatment of Pressure Injuries



Highly vapor permeable and oversized pliable quilted nylon top cover provides low shear, friction and moisture protection




Cable management system keeps cords away from the ground to help prevent trip hazards

Nylon non-slip base cover helps reduce shear and friction during bed articulation



Features & Benefits

8 adjustable comfort settings depending on the patient's therapy requirement



ALTER
STATIC

microAIR 500

- POWER
- POWER FAILURE
- LOW PRESSURE
- STATIC

ON

OFF

2-in-1 alternating pressure cycles help ensure protection for patients with low to medium risk



Soft nylon urethane backed cells provide easy cleaning and helps minimize the proliferation of bacterial growth

Cell in cell design helps prevent bottoming out in the event of a power failure

Multiple snap bladders make cells more secure with patient movement

Natural heel slope helps decrease pressure in patient's vulnerable heel area

Request a Quote or Place an Order Today!



**Contact your Inside
Sales Representative**

or



**Contact Customer Service
at 800.333.6900**

Invacare Corporation
www.invacare.com







USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

microAIR MA500 Alternating Pressure Low Air Loss Mattress System

Model #	 Dimensions	 Weight Capacity	 Flow Rate (Pump)	 Therapy Type	 Alarms	 Limited Warranty
MA500	80" L 36" W 8" H	350 lb.	5 LPM	Alternating Pressure Low Air Loss	Audio & Visual	1 Year (Cover/ Mattress) 1 Year (Power Unit)
MA500P	microAIR MA500 Pump only					
MA500M	microAIR MA500 Mattress only					

©2019

Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
®. All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise noted.
Specifications are subject
to change without
notification.
Form No. 19-022
190201

